Beam Loss Detector for new ESRF Ring

- 1) Presentation of our present BLDs with some results
- 2) What kind of BLDs & acquisition electronics to add in the future?

Questions:

- what kind of detector? (scintillator / radiator e.g. Cherenkov, Gamma?)
- what kind of visible photon to electrical convertor (PMT?)
- can all light sources (3 to 6 GeV electrons) agree on a common & optimized type of BLD ?
- this including the acquisition electronics?
- resulting in a complete & a commercial system?
- could development costs be covered by a European funding scheme?

some worries with the new L.E. ESRF Ring

- Now: no top-up, instead injections 2 or 3 times a day
 - injection efficiency reasonable
 - lifetime >50hrs (depending on filling mode)
 - degradation on (some) ID-magnets noticeable after 10 years

Future: - regular top-up,

- achievable injection efficiency ...?
- lifetime <10Hrs or even (much) lower ...?
- degradation on ID-magnets in less then a year ..?!


We think BLDs will play an important role in the future, together with dose-monitors (like the DOSFET devices)

Upgrades on the BeamLoss Detectors


The vacuum lay-out of 1 cell (1 / 32 of the Ring) with:

the 6m long ID chamber
the 2 dipoles
the 2 BeamLoss Detectors

Top View



So: in total 64 BLDs in the Ring




Beam Loss Detector:

on the inside of each of the 64 Dipoles,


Perspex Cherenkov radiator rod of 600mm x 25mm diam.


with a 10 mm tubular Lead shielding and a Photomultiplier tube for light amplification


Power & Gain Control





The PM & electronics work in integrated mode (NOT counting mode) output is filtered (1Hz BW) & amplified & digitized











Brilliance- BPM

224 Units in service for Orbit Correction



and also for Turn-by-Turn
Orbit measurements


Beam Loss Monitor

1 unit procured & installed (for 4 BLDs)




and also for Turn-by-Turn Loss measurements


on a "Fast-Compact BLD" on 1rst dipole C4




in same injection, 5 bunch mode




Injection Shot






Number of Turns



Number of Turns



Injection Shot

